



# Ctrl+Q CHEAT SHEET

Worksheet function and command reference (color code: pub, [std](#), [adv](#))

## Functions in the sheet

| Sheet name            | Description   | Extracted data                                       |
|-----------------------|---|--|
| Data                  | Original data                                       |  |
| Graph                 | Plots in binding and kinetic energies               | From data sheet                                      |
| Check                 | Plots in each scan                                  | From data sheet by "No" command                      |
| Time                  | Plots in time scale and fluence analysis            | From data sheet by "Obb" command                     |
| Fit                   | Background subtraction and peak Fitting             | From graph sheet                                     |
| Ana                   | A summary of fitting results                        | From fit sheets                                      |
| Cmp                   | A summary of BG subtracted profiles                 | From fit sheets                                      |
| Rto                   | A summary of peak areas from Ana report             | From ana sheets                                      |
| Photo                 | Plots of XAS data                                   | From data sheet                                      |
| Exp                   | Data table exported from graph sheet                | From graph sheet                                     |
| Eck                   | Data table exported from check sheet                | From check sheet                                     |
| Norm                  | Data normalized                                     | From graph sheet                                     |
| Diff                  | Data subtracted                                     | From graph sheet                                     |
| Edge                  | Pre and post-edge correction                        | From graph sheet                                     |
| Lcmb                  | Linear combination fit                              | From graph sheet                                     |
| Graph_Norm            | Plots based on the Norm sheet                       | From norm sheet                                      |
| Fit_Norm              | Peak fitting based on the Norm sheet                | From graph_norm sheet                                |
| Graph_Edge            | Plots based on the Edge sheet                       | From edge sheet                                      |
| Fit_Edge              | Peak fitting based on the Edge sheet                | From graph_edge sheet                                |
| Graph_Lcmb            | Plots based on the Lcmb sheet                       | From lcmb sheet                                      |
| Fit_Lcmb              | Peak fitting based on the Lcmb sheet                | From graph_lcmb sheet                                |
| <a href="#">Calc</a>  | <a href="#">Simulation parameters</a>               |  |
| <a href="#">Sim</a>   | <a href="#">Simulated data</a>                      | <a href="#">From Calc sheet</a>                      |
| <a href="#">Vms</a>   | <a href="#">Data table exported in VAMAS format</a> | <a href="#">From Graph sheet by "vms" or "vamas"</a> |
| <a href="#">VAMAS</a> | <a href="#">Data imported from VAMAS format</a>     |  |

## Manual data load format template

| Technique | Trigger in A1 | Queries       | Graph x-axis | Fitting         |
|-----------|---------------|---------------|--------------|-----------------|
| PES       | KE/eV         | PE & elements | BE & KE      | Yes in BE scale |
| XPS       | BE/eV         | PE & elements | BE & KE      | Yes in BE scale |
| XAS       | PE/eV         | Elements      | PE           | Yes in PE scale |

|              |       |                             |            |             |
|--------------|-------|-----------------------------|------------|-------------|
| Grating scan | GE/eV | Gap/1 <sup>st</sup> har.& e | PE         | No          |
| AES          | AE/eV | Elements                    | EE & dN/dE | No          |
| RGA          | QE/eV | NA                          | Mass       | Yes in mass |
| Manual scan  | ME/eV | NA                          | Position   | Yes in x    |
| Histogram    | HE/eV | NA                          | Position   | Yes in x    |
| Photodiode   | FE/eV | Gap/1 <sup>st</sup> har.    | PE         | No          |

### Command list in each sheet

| Command    | Cell                   | Sheet             | Result  |
|------------|------------------------|-------------------|---|
| chem       | C10                    | Graph, Cmp        | Chemical shifts   |
| elem       | C10                    | Graph             | Revise elements   |
| intp       | A1                     | Data              | Interpolate data by B1 value  |
| ana        | C10                    | Graph             | Update Fit sheet  |
| exp        | A1                     | Graph, Check, Cmp | Export data with filename-related X axis (used to paste it in the other program)    |
| exp2       | A1                     | Graph, Check, Cmp | Export data with "E/eV" labels for X axis (used to export each data into text file) |
| exp3       | A1                     | Graph             | Export data with "AE/eV" for Auger  |
| comp       | D1                     | Graph             | Compare the spectra   |
| auto       | A1                     | Graph, Cmp        | Calibrate offset and multiple factors (see the detailed functions below)            |
| cali       | A1                     | Graph             | Calibrate C1s peak position by fitting  |
| noise(n)   | A1                     | Graph             | Remove shot noise when S/N > n (default 5)  |
| ana        | D4                     | Fit               | Summarize Fit sheets  |
| ana        | A1                     | Ana               | Summarize Ana sheets into Rto sheet   |
| exp        | A1                     | Check             | Export data to be imported into elsewhere   |
| exp        | A1                     | Cmp               | Export data to be imported into elsewhere   |
|            | A1 & C1, ... = "KE/eV" | Eck               | Export text data files for each two-column  |
| debug      | A1                     | Graph             | Apply the same graph parameters to files  |
| debug      | D1                     | Fit               | Apply the same fit parameters from text files                                       |
| debugn     | A1                     | Graph             | Add reference data to be normalized   |
| debuga     | D1                     | Fit               | Apply previous fit parameters from text files                                       |
| debugf     | D1                     | Fit               | Apply RSF in the fit sheet only   |
| debugc     | D1                     | Fit               | Apply fit parameters in the fit sheet only  |
| norm, diff | A1                     | Graph             | Normalize first data by second added data   |
| cked       | A1                     | Graph             | Normalized by gold C K edge in database   |
| edge       | A1                     | Graph             | Pre and post edge correction.   |
| lcmb       | A1                     | Graph             | Linear combination of multiple spectra  |
|            |                        | Time              | Fluence analysis interpolated in any points   |
| vms        | A1                     | Graph             | Single export in vmas format  |
| vamas      | A1                     | Graph             | Multiple regions export in vmas format  |

|            |    |      |  |
|------------|----|------|--|
| phi        | A2 | Data | Export and plot from Multipak exported csv |
| multi      | A2 | Data | Plot graph from Multipak exported csv      |
| simulation | A1 |      | Simulate the spectrum with elements        |
| Imfit      | A4 | Fit  | Export Python code for Imfit               |

### Backgrounds in the fit sheet

| Type of BG                    | A1 | B1    | C1    |
|-------------------------------|----|-------|-------|
| Shirley BG                    | sh | ab/bg |       |
| Tougaard BG                   | to | ab/bg |       |
| Polynomial BG                 | po | ab/bg |       |
| Polynomial Normal BG          | po | no    | ab/bg |
| Polynomial Shirley BG         | po | sh    | ab/bg |
| Polynomial Tougaard BG        | po | to    | ab/bg |
| Polynomial Conv-Tougaard      | po | co    | ab/bg |
| Polynomial Virtual Shirley BG | po | vi    | ab/bg |
| Polynomial Edge BG            | po | ed    | ab/bg |
| Polynomial AsLS BG            | po | as    | ab/bg |
| Slope Shirley BG              | sl | sh    | ab/bg |
| Slope Tougaard BG             | sl | to    | ab/bg |
| Slope Virtual Shirley BG      | sl | vi    | ab/bg |
| Shirley Iterated BG           | sh | it    | bg    |
| Shirley Peak BG               | sh | pe    | bg    |
| Virtual Shirley BG            | vi | sh    | ab/bg |
| Tougaard Convolved            | to | co    | ab/bg |
| Arctan BG                     | ar | ab/bg |       |
| Erf BG                        | er | ab/bg |       |
| Victoreen BG                  | vi | ab/bg |       |
| Double Exponential BG         | do | ab/bg |       |
| Lognormal                     | lo | ab/bg |       |
| Sigmoid fit + spline BG       | si | fi    |       |
| Sigmoid convoluted fit        | si | co    | fi    |
| Double Sigmoid fit            | do | si    | fit   |
| User-defined function         | ud | fit   |       |
| SAXS                          | sa | fit   |       |
| CK (in fact, Arctan BG)       | ck |       |       |

### Peak shapes in the fit sheet

| Syntax             | Shape           | Option a   | Option b | #par | Ref.  |
|--------------------|-----------------|------------|----------|------|-------|
| G (0)              | Gaussian        |            |          | 3    |       |
| <i>DB G</i><br>(0) | Double Gaussian |            |          | 4    | Fityk |
| <i>EMG</i>         | Exponentially   | Distortion |          | 4    | Fityk |

|                           |                                     |                              |                            |   |                           |
|---------------------------|-------------------------------------|------------------------------|----------------------------|---|---------------------------|
|                           | Modified Gaussian                   | para.                        |                            |   |                           |
| L (1)                     | Lorentzian                          |                              |                            | 3 |                           |
| <u>DS L (1)</u>           | Doniac-Sunjic x L                   | Asymmetric para.             |                            | 5 | CasaXPS                   |
| <u>DB L (1)</u>           | Double Lorentzian                   |                              |                            | 4 | AAalyzer                  |
| <u>PEA</u>                | Pearson VII                         | Skewness                     |                            | 4 | Fityk                     |
| SGL, PGL (0-1)            | G + L, G x L (pseudo-Voigt)         |                              |                            | 5 | Unifit CasaXPS            |
| <u>ASGL, APGL</u>         | Asymmetric V Double Voigt           |                              |                            | 5 | 10.1107/S0021889884011043 |
| ESGL, EPGL                | Exponential blended Voigt           | Exponential decay parameters |                            | 5 | CasaXPS                   |
| <u>DS SGL, DS PGL</u>     | DS x L blended V                    | Asymmetric parameter         | Ratio DSL:V                | 6 | CasaXPS                   |
| <u>UG SGL, UG PGL</u>     | Ulrik Gelius blended Voigt          | Asymmetric parameter a       | Asymmetric parameter b     | 6 | CasaXPS                   |
| <u>DSV, SGL, DSV, PGL</u> | DS x Voigt blended Voigt            | Asymmetric parameter         | Ratio DSV:V                | 6 | CasaXPS                   |
| <u>TSGL</u>               | Exponential blend SGL (MultiPak)    | Tail scale                   | Tail length at half max    | 6 | MultiPak                  |
| GL (0 < shape < 1)        | G + L with the same FWHM (MultiPak) |                              |                            | 4 | MultiPak Eq. to SGL       |
| MSGL                      | 10.1002/sia.5521                    | Asymmetric parameter         | Sigmoid-center translation | 6 | 10.1002/sia.5521          |
| CGL                       | Numerical convolution G x L         |                              |                            | 4 | 10.1002/sia.2527          |
| F                         | Fano profile                        |                              |                            | 4 | 10.1103/PhysRev.124.1866  |
| FG                        | F x G                               |                              |                            | 5 |                           |
| LOGN                      | Log normal                          | Mean ( $\mu$ )               |                            | 4 |                           |

### Optimization mode of fittings

| Cell in Fit sheet              | Syntax or Font style | Optimization     |
|--------------------------------|----------------------|------------------|
| BE, FWHM, Ampl, Shape, Options | Figures with Bold    | Constraints      |
| A14                            | Solve $\chi^2$ *     | Least chi square |

|              |                         |   |
|--------------|-------------------------|---|
| A14          | Solve Abbe              | Abbe criteria   |
| A10 (EF fit) | Solve FD without Italic | Least chi square  |
| A10 (EF fit) | Solve FD with Italic    | Abbe criteria   |
| A11 (EF fit) | Solve GC without bold   | Gaussian convolution after FD + polynomial BG                             |
| A11 (EF fit) | Solve GC with bold      | FD + Polynomial BG first, Gaussian convolution together with FD + poly BG |

### Calibrations in offset/multiple factors

| A1 cell syntax in Graph sheet | Offset factor                                     | Multiple factor                          |
|-------------------------------|---|--|
| auto0                         | Set to 0  | Set to 1                                 |
| auto or auto1                 | First point to be zero                            | End point to be unity                    |
| auto10                        | Zero at point 10 from start point                 | Unity at point 10 from end point         |
| auto(1,10)                    | Zero from point 1 to 10 from start point          | Unity from point 1 to 10 from end point  |
| auto[100:101,200:201]         | Zero in BE range between 100 and 101 eV           | Unity in BE range between 200 and 201 eV |
| automax / autowf              | Zero at the lower side of a point of data         | Unity at max intensity point of data     |
| autop                         | Syntax previously done                            | Syntax previously done                   |
| auto{284.6}                   | BE at max. intensity to be calibrated in 284.6 eV | NA (BE calibration by Charging factor)   |
| auto'-7.8'                    | Charging correction at -7.8 eV for all spectra    | NA (this is based on C1s BE calibration) |
| offset10                      | Offset spectra for water fall plot                | NA                                       |

### List of element groups to be identified

| Code | Group                 | Elements to be analyzed               |
|------|-----------------------|---------------------------------------|
| AL   | Alkali metals         | Na,K,Rb,Cs                            |
| EA   | Alkaline Earth metals | Be,Mg,Ca,Sr,Ba,Ra                     |
| TM   | Transition metals     | <i>3d + 4d + 5d transition metals</i> |
| 3d   | 3d transition metals  | Sc,Ti,V,Cr,Mn,Fe,Co,Ni,Cu,Zn          |
| 4d   | 4d transition metals  | Y,Zr,Nb,Mo,Tc,Ru,Rh,Pd,Ag,Cd          |
| 5f   | 5d transition metals  | Lu,Hf,Ta,W,Re,Os,Ir,Pt,Au,Hg          |
| SM   | Semi-metals           | B,Si,Ge,As,Sb,Te                      |
| NM   | Non-metals            | C,N,O,P,S,Se                          |
| BM   | Basic metals          | Al,Ga,In,Sn,Tl,Pb,Bi                  |
| HA   | Halogens              | F,Cl,Br,I,At                          |
| NG   | Noble gases           | Ne,Ar,Kr,Xe,Rn                        |
| RM   | Rare metals           | La,Ce,Nd,Sm,Eu,Gd,Tb,Er,Tm,Yb,Th,U    |

|    |             |                               |
|----|-------------|-------------------------------|
| LA | Lanthanides | La,Ce,Nd,Sm,Eu,Gd,Tb,Er,Tm,Yb |
| AC | Actinides   | Th,U                          |

### Advanced syntax templates in the sheets

|                | Sheet | Cells | Formula         | Reference | Calibrated #1 | Calibrated #2 |
|----------------|-------|-------|-----------------|-----------|---------------|---------------|
| Extra photons  | Graph | C2    | ;100;200;333 eV |           |               |               |
| Specific scans | Graph | B8    | [1,2-4]         |           |               |               |
| BE diff        | Fit   | D14-  | (4;3)           | (4;       | 1;            | 3)            |
| Amp ratio      | Fit   | D15-  | [3.5;n3.5]      | [         | 3.5;          | n3.5]         |

Note1: "n" represents negative shift from reference.

Note2: Empty cells between brackets does not effect to the constraints.

### List of Peak area

|         | Usages  | Descriptions   | Factors to be effective   |
|---------|---|--|---|
| P. Area | Chemical state analysis   | Peak area calculated with analytical formula and without any factors   | Amplitude, FWHM   |
| S. Area | Quantification of elements under the same condition                 | Peak area normalized with atomic sensitivity factor based on photo-ionization cross-section  | Amplitude, FWHM, PE, Sensitivity based on element specified in the <i>Graph</i> sheet   |
| N. Area | Quantification of elements under the various measurement conditions | Peak area calculated in "S. Area" plus normalized with empirically calculated factors at BL CLAM2 including XPS mean-free path of photoelectrons, transmission function of electron energy analyzer based on pass energy, grating efficiency | Amplitude, FWHM, PE, KE, Sensitivity, CAE, Grating, MFP factor, a & b specified in the <i>Fit</i> sheet based on formalism from CasaXPS |

T.I./S.I./N.I. are numerically integrated areas with Trapezoidal rule applied to each corresponding area shown above.

### Batch processing for multiple files initiated by blank window

| Code | Mode             | Processing  |
|------|------------------|---|
| 1    | CLAM2 txt2xlsx   |   |
| 2    | XPS AlKa csv2txt | Multipak exported csv data to asci texts          |
| 3    | XAS SDD mca2txt  |   |
| 4    | xlsx2update      |   |
| 5    | xlsx2vamas       | Export vamas format files                         |
| 8    | PE input mode    | XPS analysis with PE input                        |
| 9    | Push charts2ppt  | All chart in each graph sheet to push them in ppt |

|    |                        |                                     |
|----|------------------------|-------------------------------------|
| 10 | CLAM2 txt2fitting      |                                     |
| 11 | Si2p fit1              | XPS fitting on Si2p low resolution  |
| 12 | Si2p fit2              | XPS fitting on Si2p high resolution |
| 13 | Au4f fit (no PE input) | XPS                                 |

### More resources

<https://github.com/heitler/xps-excel-macro>

[hideki@slri.or.th](mailto:hideki@slri.or.th)